

3-Panel Stand Instructions - Retrofitting Magnetic Stand from Tension Stand

(APPROXIMATE Overall Targets: Horizontal length 53-58", Vertical INCLUDING Legs 46-49")

Retrofit Template

You Need

You Have

Cut to

Use Instead

You Need Exact or approximate. This is basically the target length of the identified part of the stand.

You Have Based on our records, we believe you have this piece already.

Cut to In many cases, we'll tell you to cut a piece. It NORMALLY isn't so important that the cuts are precisely what we say. But it is important that the pieces are cut pretty close to identical. For example, if we say cut a 42 to 41, it would probably suffice if they were both 40.5 but not 40 and 41. If "ok", use the piece in red box without cutting.

Use Instead A substitute or additional piece. Perhaps you have 1 piece that will replace 2. Or you need to add a piece where there was none before. The support piece is not common on tension stands, for example, but is a MUST for the magnetic due to weight.

All measurements in inches

Last Modified June 2020

Additional Notes

Tension Stand Size (OLD)

27			
----	--	--	--

Upper/Lower Left

x 2

Name/Order #

27			
----	--	--	--

Upper/Lower Right

x 2

Replace tension bands with reusable zip ties.

41

Verticals

--

x 3

--

--

4

Legs

--


x 2

--

Sub Cross for T

--

Back Side View



Support

5.5

--	--	--

Parts Added

Cross	

Front to back on each leg: cap, 12", T, 12", cap



Scan QR code for latest Magnetics instructions, videos and notes. If QR does not work, enter link below into your browser.

<https://www.BrunosDraftBoards.com/Articles.asp?ID=389>



(314)961-KITS
Call us if you need help, seriously.